


# R15™ Provides Oxygen for 15 Minutes at 6L/min.

Emergency oxygen when you need it.™



## Patented Breakthrough Technology

- Provides 15 minutes of humidified oxygen at 6 liters/minute for cardiac and airway emergencies while waiting for EMS arrival
- Easy to use, no training or maintenance required
- FDA Cleared for Over The Counter; Does not require a prescription
- Requires no electric or battery power
- Unlike oxygen cylinders, R15 is non-explosive and does not require specialized storage



## The R15 can be stored safely and easily accessible in all public areas



### When would you use emergency oxygen?

- Allergic Reaction ■ Anemia ■ Anxiety/Panic Attack
- Asthma Attack ■ Carbon Monoxide Poisoning
- Cardiac Arrest ■ Concussion ■ Drowning Victim
- Drug Overdose ■ Fainting
- Shortness of Breath/Troubled Breathing


### Be prepared with emergency oxygen when you see signs of:

- Shortness of Breath
- Chest Pain
- Confusion or Slurred Speech
- Blue Lips or Nails

Call 911

### Easy to use - activates with just 3 simple steps

1




Lay unit flat. Pull tape to open mask compartment.

2



Rotate lever to line. Lever CANNOT be reversed

3



Secure mask over mouth and nose with strap.

Top 10 reasons to own the R15

Metric	Rapid Oxygen R15	Oxygen Cylinders	Portable Oxygen Concentrators
1. Non-Explosive	✓	X	✓
2. No Power Required	✓	✓	X
3. 100% Humidified Oxygen	✓	X	X
4. Emergency Use ANYWHERE (and at altitude)	✓	X	X
5. No Training/Certification Required	✓	X	✓
6. Permitted in Public Areas	✓	X	✓
7. Easy to Use	✓	X	✓
8. No Prescription Required, FDA Cleared	✓	X	X
9. Maintenance Free	✓	X	X
10. Self Contained & Portable (doesn't use ambient air)	✓	✓	X

Specifications

Flow	≥ 6 L/min
Flow Duration	≥ 15 minutes
Oxygen Volume (Minimum)	> 90 L
% Oxygen	≥ 99%
RH, Temperature	Humidified <98.6°F (37°C)
Storage Temperature Range	40°-104°F (4.4°C-40°C)
Ideal Operating Temperature Range	68°-75°F (20°-23.9°C)
Shipping Temperature Range	Not to exceed peaks above 120°F (48.8°C) or below 10°F (-12.2°C)
Shelf Life	2 years (from date of manufacture)
Weight	16 lbs (7.3Kg)
Dimensions	8" x 14" x 16" (20.32 x 35.56 x 40.64 cm)
User Interface	One Lever Operation, No Training Required
Use	Single Use - Disposable



Rapid Oxygen Company  
9 West Broad Street, Suite 550  
Stamford, CT 06902

Toll Free: (833) 727-4362  
info@rapidoxygen.com  
www.rapidoxygen.com

How to activate The R15™

